



3914 DEL AMO BLVD  
SUITE 901  
TORRANCE CA 90503  
www.celprogen.com

Phone: 310 542 8822  
Fax: 310 542 8028  
Email: info@celprogen.com  
stemcells@celprogen.com

## DATA SHEET

<b>Mouse Gastric Stem Cell pre-coated– T150 Flasks Differentiation Extracellular Matrix</b>
---

**Catalog number:** D66106-33- T150

**Description:** Mouse Gastric Stem Cell pre-coated– T150 Flasks Differentiation Extracellular Matrix pre-coated for growth, Differentiation and proliferation of Mouse Gastric Stem Cell pre-coated– T150 Flasks. The Mouse Gastric Differentiation Extracellular Matrix is available in the following cell culture ware:

**Well Plates:**

6 wells	CAT# 66106-33-6Well plate	5 per pack
12 wells	CAT# 66106-33-12Well plate	5 per pack
24 wells	CAT# 66106-33-24Well plate	5 per pack
96 wells	CAT# 66106-33-96Well plate	5 per pack
384 wells	CAT# 66106-33-384Well plate	5 per pack

**Flasks:**

T25	CAT# 66106-33-T25 flasks	10 per pack
T75	CAT# 66106-33-T75 flasks	5 per pack
T150	CAT# 66106-33-T150 flasks	5 per pack
T225	CAT# 66106-33-T225 flasks	5 per pack

**Microscope Slides:**

8 Well	CAT# 66106-33-8Well slides	5 per pack
10 Well	CAT# 66106-33-10Well slides	5 per pack
18 Well	CAT# 66106-33-18Well slides	5 per pack

**Microscope Chamber Slides:**

8 Well	CAT# 66106-33-8Well slides	5 per pack
4 Well	CAT# 66106-33-4Well slides	5 per pack
2 Well	CAT# 66106-33-2Well slides	5 per pack

**Microscope Cover Slips:**

12mm	CAT# 66106-33-CS12	10 per pack
15mm	CAT# 66106-33-CS15	10 per pack
18mm	CAT# 66106-33-CS18	10 per pack
22mm	CAT# 66106-33-CS22	10 per pack
25mm	CAT# 66106-33-CS25	10 per pack

**Petri Dishes:**

6cm	CAT# 66106-33-PD6	5 per pack
10cm	CAT# 66106-33-PD10	5 per pack



3914 DEL AMO BLVD  
SUITE 901  
TORRANCE CA 90503  
www.celprogen.com

Phone: 310 542 8822  
Fax: 310 542 8028  
Email: info@celprogen.com  
stemcells@celprogen.com

The 6wells, 12wells, 24wells, 96wells, and 384wells are available 5 well plates per pack. The T75, T150, and T225 flasks are available as 5 flasks per pack. The T25 flasks are available as 10 per pack. The 8 Well, 10 Well, and 18 Well Microscope Slides are available as 5 in pack. The 8 Well, 4 Well, and 2 Well Microscope Chamber Slides are available in 5 per pack. The 12mm, 15mm, 18mm, 22mm, and 25mm Microscope Cover Slips are available in 10 per pack. The 6cm and 10cm Petri Dishes are available in 5 per pack. For individual catalog items please visit our web site at [www.celprogen.com](http://www.celprogen.com).

**Storage Conditions:** 2 - 8°C long term storage 24 months (**thaw at room temperature prior to use**)  
Room temperature for short term storage 12 months

**Product Orders:** Before submitting an order you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

**Permits/Forms:** In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

#### **Notices**

**& Disclaimers:** *CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in humans.* The product Mouse Gastric Stem Cell pre-coated- T150 Flasks Differentiation Extracellular Matrix is established and manufactured by CELPROGEN Inc., and is for Research Use Only. The product is not for re-sale and may not be transferred to third-party prior to written request and approval by CELPROGEN Inc.